

Package ‘RDML’

January 20, 2025

Type Package

Title Importing Real-Time Thermo Cyclers (qPCR) Data from RDML Format Files

Version 1.0

LazyData true

Date 2019-06-25

Description Imports real-time thermo cyclers (qPCR) data from Real-time PCR Data Markup Language (RDML) and transforms to the appropriate formats of the 'qpcR' and 'chipPCR' packages. Contains a dendrogram visualization for the structure of RDML object and GUI for RDML editing.

License MIT + file LICENSE

URL <https://github.com/kablag/RDML>

Depends R (>= 3.2.0)

Imports checkmate (>= 1.6.2), data.table, pipeR, readxl, rlist (>= 0.4), R6 (>= 2.0.1), stringr, tools (>= 3.2), xml2 (>= 1.0), lubridate (>= 1.6.0)

Collate 'RDML.types.R' 'RDML.R' 'RDML.AsDendrogram.R' 'RDML.AsTable.R' 'RDML.GetFData.R' 'RDML.Merge.R' 'RDML.SetFData.R' 'RDML.init.R' 'functional_wrappers.R' 'rdmlEdit.R'

Suggests chipPCR, magrittr, reshape2, qpcR, dplyr, ggplot2, knitr, kfigr, MBmca, shiny, shinyjs, shinythemes, shinyMolBio, V8, testthat

RoxygenNote 6.1.1

NeedsCompilation no

Author Konstantin A. Blagodatskikh [cre, aut],
Stefan Roediger [aut],
Michal Burdukiewicz [aut] (<<https://orcid.org/0000-0001-8926-582X>>),
Andrej-Nikolai Spiess [aut]

Maintainer Konstantin A. Blagodatskikh <k.blag@yandex.ru>

Repository CRAN

Date/Publication 2019-06-25 11:40:10 UTC

Contents

adpsType	3
annotationType	4
as.character.idType	5
as.character.reactIdType	5
AsDendrogram	6
AsTable	6
baseTemperatureType	7
cdnaSynthesisMethodType	7
commercialAssayType	8
cqDetectionMethodType	9
dataCollectionSoftwareType	9
dataType	10
documentationType	11
dyeType	12
enumType	12
experimenterType	13
experimentType	13
GetFData	14
gradientType	14
idReferencesType	15
idType	16
labelFormatType	16
lidOpenType	17
loopType	17
mdpsType	18
measureType	18
MergeRDMLs	19
new	20
nucleotideType	21
oligoType	22
pauseType	22
pcrFormatType	23
primingMethodType	24
quantityType	24
quantityUnitType	25
RDML	26
RDML.AsDendrogram	29
RDML.AsTable	30
RDML.GetFData	32
RDML.SetFData	33
rdmlBaseType	34
rdmlEdit	34
rdmlIdType	35
reactIdType	35
reactType	36
runType	37

sampleType	38
sampleTypeType	39
sequencesType	40
SetFData	41
stepType	41
targetType	42
targetTypeType	43
temperatureType	43
templateQuantityType	44
thermalCyclingConditionsType	44
xRefType	45
[.GetFData	46

Index	47
--------------	-----------

adpsType	<i>adpsType R6 class.</i>
----------	---------------------------

Description

adpsType R6 class.

Usage

adpsType

Format

An [R6Class](#) generator object.

Details

Contains matrix of amplification data. Must have three columns:

eye PCR cycle at which data point was collected (every cycle must have unique number).

tmp temperature in degrees Celsius at the time of measurement (optional).

fluor raw fluorescence intensity measured.

Inherits: [rdmlBaseType](#).

Initialization

adpsType\$new(fpoints)

Fields

fpoints [assertMatrix](#). Matrix with amplification data points.

Examples

```

#cycles
cyc <- c(1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 11, 12, 13, 14, 15, 16, 17,
18, 19, 20, 21, 22, 23, 24, 25, 26, 27, 28, 29, 30, 31, 32, 33,
34, 35, 36, 37, 38, 39, 40)
#fluorescence
fluo <- c(2.0172, 2.0131, 2.0035, 2, 2.0024, 2.0056, 2.0105, 2.0179,
2.0272, 2.0488, 2.0922, 2.1925, 2.3937, 2.7499, 3.3072, 4.0966,
5.0637, 6.0621, 7.0239, 7.8457, 8.5449, 9.1282, 9.6022, 9.9995,
10.2657, 10.4989, 10.6813, 10.8209, 10.9158, 10.9668, 11.0053,
11.0318, 11.0446, 11.044, 11.0052, 10.9671, 10.9365, 10.9199,
10.897, 10.8316)
#temperature
temp <- c(55, 55, 55, 55, 54, 54, 55, 55, 55, 55, 55, 55, 55, 55, 55, 55,
55, 55, 55, 55, 55, 55, 55, 55, 56, 55, 55, 55, 55, 55, 55, 55,
55, 55, 55, 55, 55, 55, 55, 55, 55)

#combine all variables into a proper object
data <- data.frame(cyc = cyc, tmp = temp, fluor = fluo)

#create adps object
adpsType$new(data)

#create adps object without temperature data
adpsType$new(data[, -2])

```

annotationType

annotationType R6 class.

Description

Annotate samples by setting a property and its value. For example, sex could be a property with the possible values M or F. Inherits: [rdmlBaseType](#).

Usage

```
annotationType
```

Format

An [R6Class](#) generator object.

Fields

property [checkString](#). Property name

value [checkString](#). Value

Examples

```
#set sex property  
annotationType$new(property = "sex", value = "M")
```

as.character.idType *Convert idType object to character*

Description

Function to convert idType object to character.

Usage

```
## S3 method for class 'idType'  
as.character(x, ...)
```

Arguments

x	idType object.
...	Further arguments to be passed.

as.character.reactIdType
Convert reactIdType object to character

Description

Function to convert reactIdType object to character.

Usage

```
## S3 method for class 'reactIdType'  
as.character(x, ...)
```

Arguments

x	reactIdType object.
...	Further arguments to be passed.

AsDendrogram	RDML\$AsDendrogram() <i>wrapper</i>
--------------	-------------------------------------

Description

Read more at [RDML.AsDendrogram](#)

Usage

```
AsDendrogram(obj, ...)
```

Arguments

obj	RDML object.
...	AsDendrogram params.

AsTable	RDML\$AsTable() <i>wrapper</i>
---------	--------------------------------

Description

Read more at [RDML.AsTable](#)

Usage

```
AsTable(obj, ...)
```

Arguments

obj	RDML object.
...	AsTable params.

baseTemperatureType *baseTemperatureType R6 class.*

Description

Parent class for inner usage. Inherits: [rdmlBaseType](#).

Usage

baseTemperatureType

Format

An [R6Class](#) generator object.

Initialization

```
baseTemperatureType$new(duration,  
  temperatureChange = NULL, durationChange = NULL, measure = NULL, ramp =  
  NULL)
```

Fields

duration [checkCount](#). Duration of this step in seconds.

temperatureChange [checkNumber](#). Change of the temperature between two consecutive cycles:
actual temperature = temperature + (temperatureChange * cycle counter)

durationChange [checkCount](#). Change of the duration between two consecutive cycles: actual
duration = duration + (durationChange * cycle counter)

measure [measureType](#). Indicates to make a measurement and store it as meltcurve or real-time
data.

ramp [checkNumber](#). Allowed temperature change between two consecutive cycles in degrees Cel-
sius per second. If unstated, the maximal change rate is assumed.

cdnaSynthesisMethodType

cdnaSynthesisMethodType R6 class.

Description

Description of the cDNA synthesis method. Inherits: [rdmlBaseType](#).

Usage

cdnaSynthesisMethodType

Format

An [R6Class](#) generator object.

Initialization

```
cdnaSynthesisMethodType$new(enzyme = NULL,
  primingMethod = NULL, dnaseTreatment = NULL, thermalCyclingConditions =
  NULL)
```

@section Fields:

enzyme [checkString](#). Name of the enzyme used for reverse transcription.

primingMethod [primingMethodType](#).

dnaseTreatment [checkFlag](#) if TRUERNA was DNase treated prior cDNA synthesis.

thermalCyclingConditions [idReferencesType](#).

commercialAssayType *commercialAssayType R6 class.*

Description

For some commercial assays, the primer sequences may be unknown. This element allows to describe commercial assays. Inherits: [rdmlBaseType](#).

Usage

```
commercialAssayType
```

Format

An [R6Class](#) generator object.

Initialization

```
commercialAssayType$new(company, orderNumber)
```

@section Fields:

company [checkString](#).

orderNumber [checkString](#).

cqDetectionMethodType *cqDetectionMethodType R6 class.*

Description

The method used to determine the Cq value. Can take values:

"automated threshold and baseline settings"

"manual threshold and baseline settings"

"second derivative maximum"

"other"

Inherits: [enumType](#).

Usage

cqDetectionMethodType

Format

An [R6Class](#) generator object.

Initialization

cqDetectionMethodType\$new(value)

@section Fields:

value [checkString](#).

dataCollectionSoftwareType
dataCollectionSoftwareType R6 class.

Description

Software name and version used to collect and analyze the data. Inherits: [rdmlBaseType](#).

Usage

dataCollectionSoftwareType

Format

An [R6Class](#) generator object.

Initialization

```
dataCollectionSoftwareType$new(name, version)
```

@section Fields:

name [checkString](#).

version [checkString](#).

Examples

```
dataCollectionSoftwareType$new(name = "ExampleSoft",
                               version = "1.0")
```

dataType

dataType R6 class.

Description

Inherits: [rdmlBaseType](#).

Usage

```
dataType
```

Format

An [R6Class](#) generator object.

Initialization

```
dataType$new(tar, cq = NULL, excl = NULL,
             adp = NULL, mdp = NULL, endPt = NULL, bgFluor = NULL, bgFluorSlp = NULL,
             quantFluor = NULL)
```

Fields

tar [idReferencesType](#). TargetID - A reference to a target.

cq [checkNumber](#). Calculated fractional PCR cycle used for downstream quantification. Negative values express following condition: Not Available: -1.0

excl [checkString](#). Excluded. If excl is present, this entry should not be evaluated. Do not set this element to FALSE if the entry is valid. Instead, leave the entire excl element out instead. It may contain a string with a reason for the exclusion. Several reasons for exclusion should be separated by semicolons ";".

adp [adpsType](#).

mdp [mdpsType](#).

endPt [checkNumber](#). Value of the endpoint measurement.

bgFluor [checkNumber](#). Background fluorescence (the y-intercept of the baseline trend based on the estimated background fluorescence).

bgFluorSlp [checkNumber](#). Background fluorescence slope - The slope of the baseline trend based on the estimated background fluorescence. The element should be absent to indicate a slope of 0.0; If this element is present without the bgFluor element it should be ignored.

quantFluor [checkNumber](#). Quantification fluorescence - The fluorescence value corresponding to the threshold line.

Methods

```
AsDataFrame(dp.type = "adp") Represents amplification (
  dp.type = "adp"
) or melting (dp.type = "mdp") data points as data.frame
```

documentationType *documentationType R6 class.*

Description

These elements should be used if the same description applies to many samples, targets or experiments. Inherits: [rdmlBaseType](#).

Usage

```
documentationType
```

Format

An [R6Class](#) generator object.

Initialization

```
documentationType$new(id, text = NULL)
```

@section Fields:

id [idType](#). Identifier.

text [checkString](#). Text.

dyeType	<i>dyeType R6 class.</i>
---------	--------------------------

Description

Detailed information about the dye. Inherits: [rdmlBaseType](#).

Usage

```
dyeType
```

Format

An [R6Class](#) generator object.

Initialization

```
dyeType$new(id, description = NULL)
```

@section Fields:

id [idType](#). Identifier.

description [checkString](#). Description.

enumType	<i>enumType R6 class.</i>
----------	---------------------------

Description

Generic class for creating objects that can take limited list of values.

Inherits: [rdmlBaseType](#).

Usage

```
enumType
```

Format

An [R6Class](#) generator object.

Initialization

```
enumType$new(value)
```

@section Fields:

value [checkString](#). Value.

experimenterType *experimenterType R6 class.*

Description

Contact details of the experimenter. Inherits: [rdmlBaseType](#).

Usage

experimenterType

Format

An [R6Class](#) generator object.

Initialization

```
experimenterType$new(id, firstName, lastName,
  email = NULL, labName = NULL, labAddress = NULL)
```

@section Fields:

- id [idType](#). Identifier.
- firstName [checkString](#). First name.
- lastName [checkString](#). Last name.
- email [checkString](#). Email.
- labName [checkString](#). Lab name.
- labAddress [checkString](#). Lab address.

experimentType *experimentType R6 class.*

Description

A qPCR experiment. It may contain several runs ([runType](#)). Inherits: [rdmlBaseType](#).

Usage

experimentType

Format

An [R6Class](#) generator object.

Initialization

```
experimentType$new(id, description = NULL,
  documentation = NULL, run = NULL)
```

@section Fields:

id [idType](#).

description [checkString](#).

documentation list of [idReferencesType](#).

run list of [runType](#).

Methods

`AsDataFrame(dp.type = "adp", long.table = FALSE)` Represents amplification (dp.type = "adp") or melting (dp.type = "mdp") data points as data.frame. long.table = TRUE means that fluorescence data for all runs and reacts will be at one column.

GetFData	<code>RDML\$GetFData()</code> <i>wrapper</i>
----------	--

Description

Read more at [RDML.GetFData](#)

Usage

```
GetFData(obj, ...)
```

Arguments

obj	RDML object.
...	GetFData params.

gradientType	<i>gradientType R6 class.</i>
--------------	-------------------------------

Description

Details of the temperature gradient across the PCR block. Inherits: [baseTemperatureType](#).

Usage

```
gradientType
```

Format

An [R6Class](#) generator object.

Initialization

```
gradientType$new(highTemperature,  
  lowTemperature, ...)
```

Fields

highTemperature [checkNumber](#). The highest temperature of the gradient in degrees Celsius.

lowTemperature [checkNumber](#). The lowest temperature of the gradient in degrees Celsius.

... Params of parent class [baseTemperatureType](#).

<i>idReferencesType</i>	<i>idReferencesType R6 class.</i>
-------------------------	-----------------------------------

Description

Contains id of another RDML object. Inherits: [idType](#).

Usage

```
idReferencesType
```

Format

An [R6Class](#) generator object.

Initialization

```
idReferencesType$new(id)
```

Fields

id [checkString](#). Identifier.

idType	<i>idType R6 class.</i>
--------	-------------------------

Description

Contains identifier for various RDML types. Inherits: [rdmlBaseType](#).

Usage

```
idType
```

Format

An [R6Class](#) generator object.

Initialization

```
idType$new(id)
@section Fields:
id checkString. Identifier.
```

labelFormatType	<i>labelFormatType R6 class.</i>
-----------------	----------------------------------

Description

Label used for [pcrFormatType](#). Can take values:

ABC

123

A1a1

Inherits: [enumType](#).

Usage

```
labelFormatType
```

Format

An [R6Class](#) generator object.

Initialization

```
labelFormatType$new(value)
@section Fields:
value checkString.
```

lidOpenType	<i>lidOpenType R6 class.</i>
-------------	------------------------------

Description

This step waits for the user to open the lid and continues afterwards. It allows to stop the program and to wait for the user to add for example enzymes and continue the program afterwards. The temperature of the previous step is maintained. Inherits: [rdmlBaseType](#).

Usage

```
lidOpenType
```

Format

An [R6Class](#) generator object.

Initialization

```
lidOpenType$new()
```

loopType	<i>loopType R6 class.</i>
----------	---------------------------

Description

This step allows to form a loop or to exclude some steps. It allows to jump to a certain "goto" step for "repeat" times. If the "goto" step is outside of the loop range, it must have "repeat" value "0". Inherits: [rdmlBaseType](#).

Usage

```
loopType
```

Format

An [R6Class](#) generator object.

Initialization

```
loopType$new(goto, repeat.n)
```

Fields

goto [assertCount](#). The step to go to to form the loop.

repeat.n [assertCount](#). Determines how many times the loop is repeated. The first run through the loop is counted as 0, the last loop is "repeat" - 1.

mdpsType	<i>mdpsType R6 class.</i>
----------	---------------------------

Description

Contains matrix of melting data points (single data points measured during amplification).

Usage

mdpsType

Format

An [R6Class](#) generator object.

Details

Columns:

tmp (temperature in degrees Celsius at the time of measurement. Every point must have unique value.

fluor fluorescence intensity measured without any correction (including baselining).

Inherits: [rdmlBaseType](#).

Initialization

mdpsType\$new(fpoints)

@section Fields:

fpoints [assertMatrix](#). Matrix with amplification data points.

measureType	<i>measureType R6 class.</i>
-------------	------------------------------

Description

Can take values:

real time

meltcurve

Inherits: [enumType](#).

Usage

measureType

Format

An `R6Class` generator object.

Initialization

```
measureType$new(value)
```

```
@section Fields:
```

```
value checkString.
```

MergeRDMLs

Merges RDML objects

Description

Merges list of RDML objects. The first object in the list becomes base object. If experiments or runs have same name they will be combined. Reacts with same id, experiment and run overwrite each other!

Usage

```
MergeRDMLs(to.merge)
```

Arguments

to.merge RDML objects that should be merged.

Examples

```
## Not run:
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "lc96_bACTXY.rdm1", sep = "")
lc96 <- RDML$new(filename)
filename <- paste(PATH, "/extdata/", "stepone_std.rdm1", sep = "")
stepone <- RDML$new(filename)
merged <- MergeRDMLs(list(lc96,stepone))
merged$AsDendrogram()

## End(Not run)
```

new	<i>Creates new instance of RDML class object</i>
-----	--

Description

This function has been designed to import data from RDML v1.1 and v1.2 format files or from xls file generated by *Applied Biosystems 7500*. To import from xls this file have to contain Sample Setup and Multicomponent Data sheets!

Arguments

filename	string – path to file
show.progress	logical – show loading progress bar if TRUE
conditions.sep	separator for condition defined at sample name
format	string – input file format. Possible values auto, rdml, abi, excel, csv. See Details.

Details

File format options:

auto Tries to detect format by extension. .xlsx – excel, .xls – abi, .csv – csv, other – rdml

abi Reads .xls files generated by *ABI 7500 v.2*. To create such files use File>Export; check 'Sample Setup' and 'Multicomponent Data'; select 'One File'

excel .xls or .xlsx file with sheets 'description', 'adp', 'mdp'. See example file table.xlsx

csv .csv file with first column 'cyc' or 'tmp' and fluorescence data in other columns

rdml .rdml or .lc96p files

Warning

Although the format RDML claimed as data exchange format, the specific implementation of the format at devices from real manufacturers differ significantly. Currently this function is checked against RDML data from devices: *Bio-Rad CFX96*, *Roche LightCycler 96* and *Applied Biosystems StepOne*.

Author(s)

Konstantin A. Blagodatskikh <k.blag@yandex.ru>, Stefan Roediger <stefan.roediger@b-tu.de>, Michal Burdukiewicz <michalburdukiewicz@gmail.com>

Examples

```
## Not run:
## Import from RDML file
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "lc96_bACTXY.rdm1", sep = "")
lc96 <- RDML$new(filename)

## Some kind of overview for lc96
lc96$AsTable(name.pattern = sample[[react$sample$id]]$description)
lc96$AsDendrogram()

## End(Not run)
```

nucleotideType	<i>nucleotideType R6 class.</i>
----------------	---------------------------------

Description

Type of nucleic acid used as a template in the experiment. May have following values:

DNA
genomic DNA
cDNA
RNA

Usage

```
nucleotideType
```

Format

An [R6Class](#) generator object.

Details

Inherits: [enumType](#).

Initialization

```
nucleotideType$new(value)
```

@section Fields:

value [checkString](#). Value.

oligoType	<i>oligoType R6 class.</i>
-----------	----------------------------

Description

Inherits: [rdmlBaseType](#).

Usage

oligoType

Format

An [R6Class](#) generator object.

Initialization

```
oligoType$new(threePrimeTag = NULL,
              fivePrimeTag = NULL, sequence)
```

@section Fields:

threePrimeTag [checkString](#). Description of three prime modification (if present).

fivePrimeTag [checkString](#). Description of five prime modification (if present).

sequence [checkString](#).

pauseType	<i>pauseType R6 class.</i>
-----------	----------------------------

Description

This step allows to pause at a certain temperature. It is typically the last step in an amplification protocol. Inherits: [rdmlBaseType](#).

Usage

pauseType

Format

An [R6Class](#) generator object.

Initialization

```
pauseType$new(temperature)
```

Fields

temperature [checkNumber](#). The temperature in degrees Celsius maintained during the pause.

pcrFormatType *pcrFormatType R6 class.*

Description

The display format of the PCR, analogous to the the qPCR instrument run format. Inherits: [rdml-BaseType](#).

Usage

```
pcrFormatType
```

Format

An [R6Class](#) generator object.

Details

Rotor formats always have 1 column; rows correspond to the number of places in the rotor. Values for common formats are:

Format	rows	columns	rowLabel	columnLabel
single-well	1	1	123	123
48-well plate	6	8	ABC	123
96-well plate	8	12	ABC	123
384-well plate	16	24	ABC	123
1536-well plate	32	48	ABC	123
3072-well array	32	96	A1a1	A1a1
5184-well chip	72	72	ABC	123
32-well rotor	32	1	123	123
72-well rotor	72	1	123	123
100-well rotor	100	1	123	123
free format	-1	1	123	123

If rows field has value -1, the function will not try to reconstruct a plate and just display all run data in a single column. If the columns field has value 1 then the function will not display a column label.

Initialization

```
pcrFormatType$new(rows, columns, rowLabel, columnLabel)
```

@section Fields:

rows [checkCount](#).

columns [checkCount](#).

rowLabel [labelFormatType](#).

columnLabel [labelFormatType](#).

primingMethodType *primingMethodType R6 class.*

Description

The primers used in the reverse transcription. Can take values:

oligo-dt

random

target-specific

oligo-dt and random

other

Usage

primingMethodType

Format

An [R6Class](#) generator object.

Details

Inherits: [enumType](#).

Initialization

primingMethodType\$new(value)

@section Fields:

value [checkString](#). Value.

quantityType *quantityType R6 class.*

Description

A quantity is always defined by its value and its unit. Inherits: [rdmlBaseType](#).

Usage

quantityType

Format

An [R6Class](#) generator object.

Initialization

```
quantityType$new(value, unit)
@section Fields:
value checkNumber. Value.
unit quantityUnitType. Unit.
```

quantityUnitType *quantityUnitType R6 class.*

Description

The unit the quantity. Can take values:

cop copies per microliter

fold fold change

dil dilution (10 would mean 1:10 dilution)

nMol nanomol per microliter

ng nanogram per microliter

other other unit (must be linear, no exponents or logarithms allowed)

Usage

```
quantityUnitType
```

Format

An [R6Class](#) generator object.

Details

Inherits: [enumType](#).

Initialization

```
quantityUnitType$new(value)
@section Fields:
value checkString. Value.
```

RDML	<i>R6 class RDML – contains methods to read and overview fluorescence data from RDML v1.1 and v1.2 format files</i>
------	---

Description

This class is a container for RDML format data (Lefever et al. 2009). The data may be further transformed to the appropriate format of the qpcR (Ritz et al. 2008, Spiess et al. 2008) and chipPCR (Roediger et al. 2015) packages (see [RDML.new](#) for import details). Real-time PCR Data Markup Language (RDML) is the recommended file format element in the Minimum Information for Publication of Quantitative Real-Time PCR Experiments (MIQE) guidelines (Bustin et al. 2009). The inner structure of imported data faithfully reflects the structure of RDML file v1.2. All data with the exception for fluorescence values can be represented as `data.frame` by method `AsTable`. Such possibility of data representation streamlines sample filtering (by targets, types, etc.) and serves as request for `GetFData` method, which extracts fluorescence data for specified samples.

Usage

```
RDML
```

Format

An [R6Class](#) generator object.

Fields

Type, structure of data and description of fields can be viewed at RDML v1.2 file description. Names of fields are first level of XML tree.

Methods

new creates a new instance of RDML class object (see [RDML.new](#))

AsTable represent RDML data as `data.frame` (see [RDML.AsTable](#))

GetFData gets fluorescence data (see [RDML.GetFData](#))

SetFData sets fluorescence data (see [RDML.SetFData](#))

Merge merges two RDML to one (see [MergeRDMLs](#))

AsDendrogram represents structure of RDML object as dendrogram(see [RDML.AsDendrogram](#))

Author(s)

Konstantin A. Blagodatskikh <k.blag@yandex.ru>, Stefan Roediger <stefan.roediger@b-tu.de>, Michal Burdukiewicz <michalburdukiewicz@gmail.com>

References

- RDML format <http://www.rdml.org/> R6 package <http://cran.r-project.org/web/packages/R6/index.html>
- qpcR package <http://cran.r-project.org/web/packages/qpcR/index.html>
- chipPCR package: <http://cran.r-project.org/web/packages/chipPCR/index.html>
- Roediger S, Burdukiewicz M and Schierack P (2015). chipPCR: an R Package to Pre-Process Raw Data of Amplification Curves. *Bioinformatics* first published online April 24, 2015 doi:10.1093/bioinformatics/btv205
- Ritz, C., Spiess, A.-N., 2008. qpcR: an R package for sigmoidal model selection in quantitative real-time polymerase chain reaction analysis. *Bioinformatics* 24, 1549–1551. doi:10.1093/bioinformatics/btn227
- Spiess, A.-N., Feig, C., Ritz, C., 2008. Highly accurate sigmoidal fitting of real-time PCR data by introducing a parameter for asymmetry. *BMC Bioinformatics* 9, 221. doi:10.1186/1471-2105-9-221
- Bustin, S.A., Benes, V., Garson, J.A., Hellemans, J., Huggett, J., Kubista, M., Mueller, R., Nolan, T., Pfaffl, M.W., Shipley, G.L., Vandesompele, J., Wittwer, C.T., 2009. The MIQE guidelines: minimum information for publication of quantitative real-time PCR experiments. *Clin. Chem.* 55, 611–622. doi:10.1373/clinchem.2008.112797
- Lefever, S., Hellemans, J., Pattyn, F., Przybylski, D.R., Taylor, C., Geurts, R., Untergasser, A., Vandesompele, J., RDML consortium, 2009. RDML: structured language and reporting guidelines for real-time quantitative PCR data. *Nucleic Acids Res.* 37, 2065–2069. doi:10.1093/nar/gkp056

Examples

```
## EXAMPLE 1:
## internal dataset lc96_bACTXY.rdml (in 'data' directory)
## generated by Roche LightCycler 96. Contains qPCR data
## with four targets and two types.
## Import with default settings.
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "lc96_bACTXY.rdml", sep = "")
lc96 <- RDML$new(filename)

tab <- lc96$AsTable(name.pattern = paste(sample[[react$sample$id]]$description,
                                       react$id$id),
                  quantity = sample[[react$sample$id]]$quantity$value)

## Show dyes names
unique(tab$target.dyeId)
## Show types of the samples for dye 'FAM'
library(dplyr)
unique(filter(tab, target.dyeId == "FAM")$sample.type)

## Show template quantities for dye 'FAM' type 'std'#
## Not run:
COPIES <- filter(tab, target.dyeId == "FAM", sample.type == "std")$quantity
## Define calibration curves (type of the samples - 'std').
## No replicates.
library(qpcR)
CAL <- modlist(lc96$GetFData(filter(tab,
                                  target.dyeId == "FAM",
                                  sample.type == "std")),
              baseline="lin", basecyc=8:15)
```

```

## Define samples to predict (first two samples with the type - 'unkn').
PRED <- modlist(lc96$GetFData(filter(tab,
                                target.dyeId == "FAM",
                                sample.type == "unkn")),
               baseline="lin", basecyc=8:15)
## Conduct quantification.
calib(refcurve = CAL, predcurve = PRED, thresh = "cpD2",
      dil = COPIES)

## End(Not run)
## Not run:
## EXAMPLE 2:
## internal dataset lc96_bACTXY.rdml (in 'data' directory)
## generated by Roche LightCycler 96. Contains qPCR data
## with four targets and two types.
## Import with default settings.
library(chipPCR)
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "lc96_bACTXY.rdml", sep = "")
lc96 <- RDML$new(filename)

tab <- lc96$AsTable(name.pattern = paste(sample[[react$sample$id]]$description,
                                       react$id$id),
                  quantity = sample[[react$sample$id]]$quantity$value)
## Show targets names
unique(tab$target)
## Fetch cycle dependent fluorescence for HEX chanel
tmp <- lc96$GetFData(filter(tab, target == "bACT", sample.type == "std"))
## Fetch vector of dillutions
dilution <- filter(tab, target.dyeId == "FAM", sample.type == "std")$quantity

## Use plotCurves function from the chipPCR package to
## get an overview of the amplification curves
tmp <- as.data.frame(tmp)
plotCurves(tmp[,1], tmp[,-1])
par(mfrow = c(1,1))
## Use inder function from the chipPCR package to
## calculate the Cq (second derivative maximum, SDM)
SDMout <- sapply(2L:ncol(tmp), function(i) {
  SDM <- summary(inder(tmp[, 1], tmp[, i]), print = FALSE)[2]
})

## Use the effcalc function from the chipPCR package and
## plot the results for the calculation of the amplification
## efficiency analysis.
plot(effcalc(dilution, SDMout), CI = TRUE)

## End(Not run)
## Not run:
## EXAMPLE 3:
## internal dataset BioRad_qPCR_melt.rdml (in 'data' directory)
## generated by Bio-Rad CFX96. Contains qPCR and melting data.
## Import with custom name pattern.

```

```
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "BioRad_qPCR_melt.rdml", sep = "")
cfx96 <- RDML$new(filename)
## Use plotCurves function from the chipPCR package to
## get an overview of the amplification curves
library(chipPCR)
## Extract all qPCR data
tab <- cfx96$AsTable()
cfx96.qPCR <- as.data.frame(cfx96$GetFData(tab))
plotCurves(cfx96.qPCR[,1], cfx96.qPCR[,-1], type = "l")

## Extract all melting data
cfx96.melt <- cfx96$GetFData(tab, dp.type = "mdp")
## Show some generated names for samples.
colnames(cfx96.melt)[2L:5]
## Select columns that contain
## samples with dye 'EvaGreen' and have type 'pos'
## using filtering by names.
cols <- cfx96$GetFData(filter(tab, grepl("pos_EvaGreen$", fdata.name)),
                        dp.type = "mdp")
## Conduct melting curve analysis.
library(qpcR)
invisible(meltcurve(cols, fluos = 2:ncol(cols),
                    temps = rep(1, ncol(cols) - 1)))

## End(Not run)
```

RDML.AsDendrogram

Represents structure of RDML file as dendrogram

Description

Plots and/or returns the structure of RDML file as [dendrogram](#) (tree-like structure.)

Arguments

plot.dendrogram
plots dendrogram if TRUE

Value

dendrogram object

Author(s)

Konstantin A. Blagodatskikh <k.blag@yandex.ru>, Stefan Roediger <stefan.roediger@b-tu.de>, Michal Burdukiewicz <michalburdukiewicz@gmail.com>

Examples

```
## Not run:
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "BioRad_qPCR_melt.rdml", sep = "")
cfx96 <- RDML$new(filename)
#plot dendrogram
cfx96$AsDendrogram()
#assign dendrogram to the object
dendr <- cfx96$AsDendrogram(plot.dendrogram = FALSE)

## End(Not run)
```

RDML.AsTable

Represents fields of RDML object as data.frame

Description

Formats particular fields of RDML object as data.frames, filters or passes them to [RDML.GetFData](#) and [RDML.SetFData](#) functions.

Arguments

<code>.default</code>	list of default columns
<code>name.pattern</code>	expression to form fdata.name (see Examples)
<code>add.columns</code>	list of additional columns
<code>treat.null.as.na</code>	if value is NULL then convert it to NA. Helps to deal with incomplete records.
<code>...</code>	additional columns

Details

By default input this function forms data.frame with following columns:

exp.id experiment\$*id*
run.id run\$*id*
react.id react\$*id*
position react\$*position*
sample react\$*sample*
target data\$*target*
target.dyeId target[[data\$*target*]]\$*dyeId*
sample.type sample[[react\$*sample*]]\$*type*

You can overload default columns list by parameter `.default` but note that columns

`exp.id`, `run.id`, `react.id`, `target`

are necessary for usage AsTable output as input for GetFData and SetFData. Additional columns can be introduced by specifying them at input parameter . . . (see Examples). All default and additional columns accession expressions must be named.

Experiment, run, react and data to which belongs each fluorescence data vector can be accessed by experiment, run, react, data (see Examples).

Result table does not contain data from experiments with ids starting with '!'.

Author(s)

Konstantin A. Blagodatskikh <k.blag@yandex.ru>, Stefan Roediger <stefan.roediger@b-tu.de>, Michal Burdukiewicz <michalburdukiewicz@gmail.com>

Examples

```
## Not run:
## internal dataset stepone_std.rdml (in 'data' directory)
## generated by Applied Biosystems Step-One. Contains qPCR data.
library(chipPCR)
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "stepone_std.rdml", sep = "")
stepone <- RDML$new(filename)
## Mark fluorescence data which Cq > 30 and add quantities to
## AsTable output.
## Names for fluorescence data will contain sample name and react
## positions
tab <- stepone$AsTable(
  name.pattern = paste(react$sample$id, react$position),
  add.columns = list(cq30 = if(data$cq >= 30) ">=30" else "<30",
    quantity = sample[[react$sample$id]]$quantity$value)
)
## Show cq30 and quantities
tab[, c("cq30", "quantity")]
## Get fluorescence values for 'std' type samples
## in format ready for ggplot function
library(dplyr)
fdata <- stepone$GetFData(
  filter(tab, sample.type == "std"),
  long.table = TRUE)
## Plot fdata with colour by cq30 and shape by quantity
library(ggplot2)
ggplot(fdata, aes(x = cyc, y = fluor,
  group = fdata.name,
  colour = cq30,
  shape = as.factor(quantity))) +
  geom_line() + geom_point()

## End(Not run)
```

`RDML.GetFData`*Gets fluorescence data vectors from RDML object*

Description

Gets fluorescence data vectors from RDML object for specified method of experiment.

Arguments

<code>request</code>	Output from <code>AsTable</code> method(RDML.AsTable)
<code>dp.type</code>	Type of fluorescence data (i.e. 'adp' for qPCR or 'mdp' for melting)
<code>long.table</code>	Output table is ready for ggplot (See RDML.AsTable for example)

Value

matrix which contains selected fluorescence data and additional information from request if `long.table = TRUE`.

Author(s)

Konstantin A. Blagodatskikh <k.blag@yandex.ru>, Stefan Roediger <stefan.roediger@b-tu.de>, Michal Burdukiewicz <michalburdukiewicz@gmail.com>

Examples

```
## Not run:
## internal dataset BioRad_qPCR_melt.rdml (in 'data' directory)
## generated by Bio-Rad CFX96. Contains qPCR and melting data.
## Import without splitting by targets/types and with
## custom name pattern.
PATH <- path.package("RDML")
filename <- paste(PATH, "/extdata/", "BioRad_qPCR_melt.rdml", sep = "")
cfx96 <- RDML$new(filename)
## Select melting fluorescence data with sample.type 'unkn'.
library(dplyr)
tab <- cfx96$AsTable()
fdata <- cfx96$GetFData(filter(tab, sample.type == "unkn"),
                        dp.type = "adp")
## Show names for obtained fdata
colnames(fdata)

## End(Not run)
```

RDML.SetFData	<i>Sets fluorescence data vectors to RDML object</i>
---------------	--

Description

Sets fluorescence data vectors to RDML object for specified method of experiment.

Arguments

data	matrix containing in the first column data corresponding to all fluorescence values in the following columns. The name of the first column is the name of variable and names of other column are <code>fdata.names</code> (links to rows at description).
description	output from <code>AsTable</code> function that describes fluorescence data.
fdata.type	'adp' for qPCR, 'mdp' for melting data.

Examples

```
## Not run:
PATH <- path.package("RDML")
filename <- paste0(PATH, "/extdata/", "stepone_std.rdml")
cfx96 <- RDML$new(filename)
## Use plotCurves function from the chipPCR package to
## get an overview of the amplification curves
library(chipPCR)
## Extract all qPCR data
tab <- cfx96$AsTable()
tab2 <- tab
tab2$run.id <- "cpp"
cfx96.qPCR <- as.data.frame(cfx96$GetFData(tab))
cpp <- cbind(cyc = cfx96.qPCR[, 1],
  apply(cfx96.qPCR[, -1], 2,
    function(y) CPP(x = cfx96.qPCR[, 1], y = y)$y.norm))
cfx96$SetFData(cpp, tab2)
library(ggplot2)
library(gridExtra)
cfx96.gg <- cfx96$GetFData(tab, long.table = TRUE)
cpp.gg <- cfx96$GetFData(tab2,
  long.table = TRUE)
plot1 <- ggplot(cfx96.gg, aes(x = cyc, y = fluor,
  group=fdata.name)) +
  geom_line() +
  ggtitle("Raw data")
plot2 <- ggplot(cpp.gg, aes(x = cyc, y = fluor,
  group=fdata.name)) +
  geom_line() +
  ggtitle("CPP processed data")
grid.arrange(plot1, plot2, nrow=2)

## End(Not run)
```

rdmlBaseType	<i>Base R6 class for RDML package.</i>
--------------	--

Description

Most classes from RDML package inherit this class. It is designed for internal usage and should not be directly accessed.

Usage

```
rdmlBaseType
```

Format

An [R6Class](#) generator object.

Initialization

```
rdmlBaseType$new()
```

Methods

`.asXMLnodes(node.name)` Represents object as XML nodes. Should not be called directly. `node.name` – name of the root node for the generated XML tree.

`print(...)` prints object

rdmlEdit	<i>RDML Editor Graphical User Interface</i>
----------	---

Description

Launches graphical user interface that can edit RDML metadata and show qPCR or melting curves.

Usage

```
rdmlEdit()
```

rdmlIdType	<i>rdmlIdType R6 class.</i>
------------	-----------------------------

Description

This element can be used to assign a publisher and id to the RDML file.
Inherits: [rdmlBaseType](#).

Usage

```
rdmlIdType
```

Format

An [R6Class](#) generator object.

Initialization

```
rdmlIdType$new(publisher, serialNumber,  
MD5Hash = NULL)
```

Fields

publisher [checkString](#). RDML file publisher.

serialNumber [checkString](#). Serial number.

MD5Hash [checkString](#). An MD5Hash calculated over the complete file after removing all rdmlID-Types and all whitespaces between elements.

reactIdType	<i>reactIdType R6 class.</i>
-------------	------------------------------

Description

Contains identifier for reactType. Inherits: [rdmlBaseType](#).

Usage

```
reactIdType
```

Format

An [R6Class](#) generator object.

Initialization

```
reactIdType$new(id)
@section Fields:
id checkCount. Identifier.
```

reactType	<i>reactType R6 class.</i>
-----------	----------------------------

Description

A reaction is an independent chemical reaction corresponding for example to a well in a 96 well plate, a capillary in a rotor, a through-hole on an array, etc. Inherits: [rdmlBaseType](#).

Usage

```
reactType
```

Format

An [R6Class](#) generator object.

Details

The ID of this reaction

Schemas :

- rotor : assign IDs according to the position of the sample on the rotor (1 for the 1st sample, 2 for the 2nd, ...)
- plate (96/384/1536 well) : the IDs are assigned in a row-first/column-second manner. For each row, the samples are numbered according to the increasing column number. At the end of a row, the numbering starts at the first column of the next row. An example for this type of plate can be found below :

	1	2	3	...
A	1	2	3	
B	13	14		
...				

or

	1	2	3	...
1	1	2	3	
2	13	14		
...				

- multi-array plate (BioTrove) : the IDs are assigned in a row-first/column-second manner, ignoring the organisation of sub-arrays. For each row, the samples are numbered according to the increasing column number. At the end of a row, the the next row. An example for this type of plate can be found below : todo...

Initialization

```
reactType$new(id, sample, data = NULL, pcrFormat = pcrFormatType$new(8, 12, labelFormatType$new("123"))
```

@section Fields:

id [reactIdType](#). See 'Details'.

sample [idReferencesType](#). SampleID - A reference to a sample.

data list of [dataType](#).

position Human readable form of the react id (i.e. '13' -> 'B1')..

Methods

`AsDataFrame(dp.type = "adp")` Represents amplification (dp.type = "adp") or melting (dp.type = "mdp") data points of all targets as one data.frame

`.recalcPosition(pcrformat)` Converts react id to the human readable form (i.e. '13' -> 'B1').

This converted value can be accessed by position field. pcrFormat is pcrFormatType. Currently, only 'ABC' and '123' are supported as labels. For '123' '123' the Position will look like 'r01c01', for 'ABC' '123' it will be 'A01' and for '123' 'ABC' it will be 01A. 'ABC' 'ABC' is not currently supported. Note that 'ABC' will result in loss of information if the experiment contains more than 26 rows!

runType

runType R6 class.

Description

A run is a set of reactions performed in one "run", for example one plate, one rotor, one array, one chip. Inherits: [rdmlBaseType](#).

Usage

```
runType
```

Format

An [R6Class](#) generator object.

Initialization

```
runType$new(id, description = NULL,
  documentation = NULL, experimenter = NULL, instrument = NULL,
  dataCollectionSoftware = NULL, backgroundDeterminationMethod = NULL,
  cqDetectionMethod = NULL, thermalCyclingConditions = NULL, pcrFormat,
  runDate = NULL, react = NULL)
```

Fields

id [idType](#).
 description [checkString](#).
 documentation list of [idReferencesType](#).
 experimenter list of [idReferencesType](#).
 instrument [checkString](#). Description of the instrument used to acquire the data.
 dataCollectionSoftware [dataCollectionSoftwareType](#). Description of the software used to analyze/collect the data.
 backgroundDeterminationMethod [checkString](#). Description of method used to determine the background.
 cqDetectionMethod [cqDetectionMethodType](#). Description of method used to calculate the quantification cycle.
 thermalCyclingConditions [idReferencesType](#). The program used to acquire the data.
 pcrFormat [adpsType](#).
 runDate [adpsType](#). Time stamp of data acquisition.
 react list of [adpsType](#).

Methods

AsDataFrame(dp.type = "adp") Represents amplification (dp.type = "adp") or melting (dp.type = "mdp") data points as data.frame

sampleType

sampleType R6 class.

Description

A sample is a template solution with defined concentration. Since dilutions of the same material differ in concentration, they are considered different samples. A technical replicate samples should contain the same name (reactions are performed on the same material), and biological replicates should contain different names (the template derived from the different biological replicates is divergent). Serial dilutions in a standard curve must have different names (preferably stating their dilution). Inherits: [rdmlBaseType](#).

Usage

sampleType

Format

An [R6Class](#) generator object.

Initialization

```
sampleType$new(id, description = NULL,
  documentation = NULL, xRef = NULL, annotation = NULL, type =
  sampleTypeType$new("unkn"), interRunCalibrator = FALSE, quantity = NULL,
  calibratorSample = FALSE, cdnaSynthesisMethod = NULL, templateQuantity =
  NULL)
```

@section Fields:

id [idType](#). Concentration of the template in nanogram per microliter in the final reaction mix.

description [checkString](#).

documentation list of [idReferencesType](#).

xRef list of [xRefType](#).

annotation list of [annotationType](#).

type [sampleTypeType](#).

interRunCalibrator [checkFlag](#). TRUE if this sample is used as inter run calibrator.

quantity [quantityType](#). Quantity - The reference quantity of this sample. It should be only used if the sample is part of a standard curve. The provided value will be used to quantify unknown samples in absolute quantification assays. Only the use of positive integers (like 1, 10, 100, 1000) and fractions (e.g. 1, 0.1, 0.01, 0.001) is acceptable. The use of exponents (1, 2, 3, 4 or -1, -2, -3, -4) is forbidden, because it will not be interpreted as 10E1, 10E2, 10E3, 10E4 or 10E-1, 10E-2, 10E-3, 10E-4.

calibratorSample [checkFlag](#). TRUE if this sample is used as calibrator sample.

cdnaSynthesisMethod [cdnaSynthesisMethodType](#).

templateQuantity [templateQuantityType](#).

sampleTypeType

sampleTypeType R6 class.

Description

Can take values:

unkn unknown sample

ntc non template control

nac no amplification control

std standard sample

ntp no target present

nrt minusRT

pos positive control

opt optical calibrator sample

Usage

sampleTypeType

Format

An [R6Class](#) generator object.

Details

Inherits: [enumType](#).

Initialization

```
sampleTypeType$new(value)
```

@section Fields:

value [checkString](#). Value.

sequencesType	<i>sequencesType R6 class.</i>
---------------	--------------------------------

Description

Inherits: [rdmlBaseType](#).

Usage

sequencesType

Format

An [R6Class](#) generator object.

Initialization

```
sequencesType$new(forwardPrimer = NULL,  
reversePrimer = NULL, probe1 = NULL, probe2 = NULL, amplicon = NULL)
```

@section Fields:

forwardPrimer [oligoType](#).

reversePrimer [oligoType](#).

probe1 [oligoType](#).

probe2 [oligoType](#).

amplicon [oligoType](#).

SetFData	RDML\$SetFData() <i>wrapper</i>
----------	---------------------------------

Description

Read more at [RDML.SetFData](#)

Usage

```
SetFData(obj, ...)
```

Arguments

obj	RDML object.
...	SetFData params.

stepType	<i>stepType R6 class.</i>
----------	---------------------------

Description

Inherits: [rdmlBaseType](#).

Usage

```
stepType
```

Format

An [R6Class](#) generator object.

Initialization

```
stepType$new(nr, description = NULL,
  temperature = NULL, gradient = NULL, loop = NULL, pause = NULL, lidOpen =
  NULL)
```

Fields

nr [checkCount](#). The incremental number of the step. First step should have value 1. The increment between steps should be constant and equivalent to 1.

description [checkString](#).

temperature [temperatureType](#).

gradient [gradientType](#).

loop [loopType](#).
 pause [pauseType](#).
 lidOpen [lidOpenType](#).

targetType	<i>targetType R6 class.</i>
------------	-----------------------------

Description

A target is a PCR reaction with defined set of primers. PCR reactions for the same gene with distinct primer sequences are considered different targets. Inherits: [rdmlBaseType](#).

Usage

targetType

Format

An [R6Class](#) generator object.

Initialization

```
targetType$new(id, description = NULL,
  documentation = NULL, xRef = NULL, type, amplificationEfficiencyMethod =
  NULL, amplificationEfficiency = NULL, amplificationEfficiencySE = NULL,
  detectionLimit = NULL, dyeId, sequences = NULL, commercialAssay = NULL)
```

Fields

id [idType](#).
 description [checkString](#).
 documentation list of [idReferencesType](#).
 xRef list of [xRefType](#).
 type [targetTypeType](#).
 amplificationEfficiencyMethod [checkString](#).
 amplificationEfficiency [checkNumber](#).
 amplificationEfficiencySE [checkNumber](#).
 detectionLimit [checkNumber](#).
 dyeId [idReferencesType](#).
 sequences [sequencesType](#).
 commercialAssay [commercialAssayType](#).

targetTypeType	<i>targetTypeType R6 class.</i>
----------------	---------------------------------

Description

Can take values:

ref reference target

toi target of interest

Inherits: [enumType](#).

Usage

targetTypeType

Format

An [R6Class](#) generator object.

Initialization

targetTypeType\$new(value)

@section Fields:

value [checkString](#).

temperatureType	<i>temperatureType R6 class.</i>
-----------------	----------------------------------

Description

This step keeps a constant temperature on the heat block. Inherits: [baseTemperatureType](#).

Usage

temperatureType

Format

An [R6Class](#) generator object.

Initialization

temperatureType\$new(temperature, ...)

Fields

temperature [checkNumber](#). The temperature of the step in degrees Celsius.
 ... Params of parent class [baseTemperatureType](#).

templateQuantityType *templateQuantityType R6 class.*

Description

Inherits: [rdmlBaseType](#).

Usage

templateQuantityType

Format

An [R6Class](#) generator object.

Initialization

templateQuantityType\$new(conc, nucleotide)

@section Fields:

conc [checkNumber](#). Concentration of the template in nanogram per microliter in the final reaction mix.

nucleotide [nucleotideType](#).

thermalCyclingConditionsType
thermalCyclingConditionsType R6 class.

Description

A cycling program for PCR or to amplify cDNA. Inherits: [rdmlBaseType](#).

Usage

thermalCyclingConditionsType

Format

An [R6Class](#) generator object.

Initialization

```
thermalCyclingConditionsType$new(id,
  description = NULL, documentation = NULL, lidTemperature = NULL,
  experimenter = NULL, step)
```

Fields

id [idType](#).

description [checkString](#).

documentation list of [idReferencesType](#).

lidTemperature [checkNumber](#). The temperature in degrees Celsius of the lid during cycling.

experimenter list of [idReferencesType](#). Reference to the person who made or uses this protocol.

step list of [stepType](#). The steps a protocol runs through to amplify DNA.

xRefType

xRefType R6 class.

Description

Inherits: [rdmlBaseType](#).

Usage

```
xRefType
```

Format

An [R6Class](#) generator object.

Initialization

```
xRefType$new(name = NULL, id = NULL)
```

@section Fields:

name [checkString](#). Reference to an external database, for example "GenBank".

id [checkString](#). The ID of the entry within the external database, for example "AJ832138".

`[.GetFData`*Extract data points from RDML object*

Description

Extract data points from RDML object as.data.frame.

Usage

```
## S3 method for class 'RDML'  
x[i, j, dp.type = "adp"]
```

Arguments

<code>x</code>	RDML object.
<code>i, j</code>	indices.
<code>dp.type</code>	Type of fluorescence data (i.e. 'adp' for qPCR or 'mdp' for melting).

Index

- * **Bio-Rad**
 - RDML, 26
- * **CFX96**
 - RDML, 26
- * **IO**
 - RDML, 26
- * **LightCycler**
 - RDML, 26
- * **RDML**
 - RDML, 26
- * **StepOne**
 - RDML, 26
- * **datasets**
 - adpsType, 3
 - annotationType, 4
 - baseTemperatureType, 7
 - cdnaSynthesisMethodType, 7
 - commercialAssayType, 8
 - cqDetectionMethodType, 9
 - dataCollectionSoftwareType, 9
 - dataType, 10
 - documentationType, 11
 - dyeType, 12
 - enumType, 12
 - experimenterType, 13
 - experimentType, 13
 - gradientType, 14
 - idReferencesType, 15
 - idType, 16
 - labelFormatType, 16
 - lidOpenType, 17
 - loopType, 17
 - mdpsType, 18
 - measureType, 18
 - nucleotideType, 21
 - oligoType, 22
 - pauseType, 22
 - pcrFormatType, 23
 - primingMethodType, 24
 - quantityType, 24
 - quantityUnitType, 25
 - rdmlBaseType, 34
 - rdmlIdType, 35
 - reactIdType, 35
 - reactType, 36
 - runType, 37
 - sampleType, 38
 - sampleTypeType, 39
 - sequencesType, 40
 - stepType, 41
 - targetType, 42
 - targetTypeType, 43
 - temperatureType, 43
 - templateQuantityType, 44
 - thermalCyclingConditionsType, 44
 - xRefType, 45
- * **file**
 - RDML, 26
- * **hplot**
 - rdmlEdit, 34
- * **manip**
 - [.GetFData, 46
 - as.character.idType, 5
 - as.character.reactIdType, 5
 - RDML.AsDendrogram, 29
 - RDML.AsTable, 30
 - RDML.GetFData, 32
- * **qPCR**
 - RDML, 26
 - [.GetFData, 46
 - [.RDML ([.GetFData), 46
 - adpsType, 3, 10, 38
 - annotationType, 4, 39
 - as.character.idType, 5
 - as.character.reactIdType, 5
 - AsDendrogram, 6
 - assertCount, 17
 - assertMatrix, 3, 18

- AsTable, 6
- baseTemperatureType, 7, 14, 15, 43, 44
- cdnaSynthesisMethodType, 7, 39
- checkCount, 7, 23, 36, 41
- checkFlag, 8, 39
- checkNumber, 7, 10, 11, 15, 22, 25, 42, 44, 45
- checkString, 4, 8–16, 19, 21, 22, 24, 25, 35, 38–43, 45
- commercialAssayType, 8, 42
- cqDetectionMethodType, 9, 38
- dataCollectionSoftwareType, 9, 38
- dataType, 10, 37
- dendrogram, 29
- documentationType, 11
- dyeType, 12
- enumType, 9, 12, 16, 18, 21, 24, 25, 40, 43
- experimenterType, 13
- experimentType, 13
- GetFData, 14
- gradientType, 14, 41
- idReferencesType, 8, 10, 14, 15, 37–39, 42, 45
- idType, 11–15, 16, 38, 39, 42, 45
- labelFormatType, 16, 23
- lidOpenType, 17, 42
- loopType, 17, 42
- mdpsType, 10, 18
- measureType, 7, 18
- MergeRDMLs, 19, 26
- new, 20
- nucleotideType, 21, 44
- oligoType, 22, 40
- pauseType, 22, 42
- pcrFormatType, 16, 23
- primingMethodType, 8, 24
- quantityType, 24, 39
- quantityUnitType, 25, 25
- R6Class, 3, 4, 7–13, 15–19, 21–26, 34–38, 40–45
- RDML, 26
- RDML.AsDendrogram, 6, 26, 29
- RDML.AsTable, 6, 26, 30, 32
- RDML.GetFData, 14, 26, 30, 32
- RDML.new, 26
- RDML.new (new), 20
- RDML.SetFData, 26, 33, 41
- rdmlBaseType, 3, 4, 7–13, 16–18, 22–24, 34, 35–38, 40–42, 44, 45
- rdmlEdit, 34
- rdmlIdType, 35
- reactIdType, 35, 37
- reactType, 36
- runType, 13, 14, 37
- sampleType, 38
- sampleTypeType, 39, 39
- sequencesType, 40, 42
- SetFData, 41
- stepType, 41, 45
- targetType, 42
- targetTypeType, 42, 43
- temperatureType, 41, 43
- templateQuantityType, 39, 44
- thermalCyclingConditionsType, 44
- xRefType, 39, 42, 45