

# The *IdMappingRetrieval* package in Bioconductor: Collecting and caching identifier mappings from online sources.

Alex Lisovich <sup>†</sup>, Roger S. Day <sup>†‡</sup>

April 16, 2015

<sup>†</sup>Department of Biomedical Informatics, <sup>†‡</sup> Department of Biostatistics,  
University of Pittsburgh

## 1 Overview

Research which integrates data from multiple data platforms must of course merge on samples processed in parallel on the platforms. However, exploiting the full biological significance of the data depends on merging on the respective features as well. The features have platform-specific biological identifiers, so identifier mapping is critical to this merging. The *IdMappingRetrieval* package allows initial acquisition of biological identifier mappings (ID maps) from online bioinformatics services, with caching in localdata repositories for subsequent fast retrieval. An ID map is a one-to-many map from one ID type (called the primary key) to another (called the secondary key). The services currently supported are NetAffx, DAVID, and Ensembl. The package employs a unified interface for accessing these services, so that the local repositories will be easy to create, update, and use. Aside from identifier maps themselves, the the service's complete annotation data sets are also accessible through the same mechanism. . Therefore, although ID mapping is the primary goal, secondarily the package performs as a generic annotation data collection tool if desired.

The objects produced by this package are specifically suited for use by the package *IdMappingAnalysis*, currently in preparation for Bioconductor. The purpose of *IdMappingAnalysis* is to characterize and compare two or more ID maps.

Retrieval is typically subsetting by Affymetrix array.

## 2 The package architecture

### 2.1 Two types of services: query-based and file-based

The *IdMappingRetrieval* package retrieves data using service objects. The package supports two types of services, query-based and file-based. The query-based services support automated retrieval either through a standard service API, through third-party R packages, or through simulated interaction with service web pages. The file-based services assist the manual download of raw data, followed by extraction of ID maps and other information. For diverse reasons they generally require more active user participation.

#### 1. Query-based services

- (a) Affymetrix annotation file repository [1],  
Data source: <http://www.bioconductor.org/packages/2.2/bioc/html/AffyCompatible.html>.  
Service class name: `AnnotationAffx`.  
Access method: *AffyCompatible* Bioconductor package.

- (b) DAVID online query system [2,3],  
Data source: <http://www.bioconductor.org/packages/release/bioc/html/DAVIDQuery.html>  
Service class name: AnnotationDavid.  
Access method: *DAVIDQuery* Bioconductor package.
- (c) Ensembl online query system [4,5], accessed using biomaRt R package  
Data source: <http://www.bioconductor.org/packages/2.2/bioc/html/biomaRt.html>  
Service class name: AnnotationEnsembl.  
Access method: *biomaRt* Bioconductor package.
- (d) EnVision online query system [6], accessed using the SOAP-based Web Services.  
Data source: **LINK NEEDED**.  
Service class name: AnnotationEnvision  
Access method: *ENVISIONQuery* Bioconductor package.

## 2. File-based services

- (a) Affymetrix NetAffx batch system supporting data download in chunks [1].  
Data source: <http://www.affymetrix.com/analysis/index.affx>.  
Service class name: AnnotationNetxAffx.  
Access method: See Appendix A.
- (b) DAVID knowledge base data repository [2,3].  
Data source: [http://david.abcc.ncifcrf.gov/knowledgebase/DAVID\\_knowledgebase.html](http://david.abcc.ncifcrf.gov/knowledgebase/DAVID_knowledgebase.html).  
Service class name: AnnotationDavidCsv.  
Access method: See Appendix B.
- (c) Ensembl download service [5,6].  
Data source: <http://uswest.ensembl.org>.  
Service class name: AnnotationEnsemblCsv..  
Access method: See Appendix C.

## 2.2 User tasks, services, classes and objects

The package utilizes the R.oo framework for S3 class based development ([7],

<http://cran.r-project.org/web/packages/R.oo/index.html>), which we found to be well suited for the task. As seen in the table above, different subclasses of the Annotation class support access to each specific ID map retrieval service.

1. For each service of interest:
  - (a) If necessary, set user credentials for login.
  - (b) Instantiate the service object as an instance of the Annotations subclass suitable for the service of interest, Specify the query parameters.
  - (c) Invoke a data collection method (either getIdList or getDataFrame) on the service object. or a set of objects each of which represents the particular service..
2. For a list of services:
  - (a) Loop over the services as above.
  - (b) Alternatively, for convenience, instantiate the ServiceManager and use corresponding methods.

Some parameters are common for all service object types.

- `cacheFolderName` - a subdirectory name where data for given service object will be cached (“EnVision”, “DAVID” etc.)
- `primaryColumn` and `secondaryColumn` - the column names which a service object uses to extract the Id Map from the raw data frame. In case `primaryColumn`/`secondaryColumn` contain multiple column names, the service object merges the corresponding raw data frame columns into one (service specific)
- `swap` - a boolean. If true, the `primary`/`secondary` columns will be swapped at the end of the data retrieval process. Swapping reverses the ‘one to many’ identifier relation
- `species` - the species the collected data should be subsetted on

The user can override default values during the service object’s creation. The default values are typically suitable. For the file-based services, the constructor has an additional parameter:

- the location of the file(s) or the directory containing the input data.

Note that micro array type is not an attribute of a service object and is supplied as an argument to the data retrieval method, which allows to reuse the same service object for multiple micro array types data collection. It also allows to define micro array type interactively in case array type argument is equal to ‘menu’ (see help pages for details).

Class	Data Source	Query versus File-based
AnnotationAffx (based on AffyCompatible)	Affymetrix annotation file repository	Query
AnnotationNetAffx	Affymetrix NetAffx batch query system	File (Appendix A)
AnnotationDavid (based on DAVIDQuery)	DAVID online query system	Query
AnnotationDavidCsv	DAVID knowledge base	File (Appendix B)
AnnotationEnsembl (based on biomaRt)	Ensembl	Query
AnnotationEnsemblCsv	Ensembl front end	File (Appendix C)
AnnotationEnvision (based on ENVISIONQuery)	EnVision Web Services	Query

## 2.3 The caching subsystem

Due to the large amount of data available for high density micro arrays, the process of online data retrieval takes significant time ranging from a few minutes for a “fast” service to up to an hour (as in case of EnVision Web Services). It would be impractical to retrieve the same data again and again at the beginning of each analysis session given the fact that the data in online repository typically get updated few times a year. The IdMappingRetrieval caching mechanism allows the system to store the retrieved data locally allowing to reuse them in subsequent processing sessions significantly (hundreds of times) increasing the retrieval speed. When request for data retrieval is received, the system checks if the data requested already stored in caching directory and initiates the online retrieval process only when data of interest is not available. The system can be forced to retrieve the data online even if present in caching subsystem thus allowing to update the cached data if necessary. The caching subsystem is implemented as a hierarchical set of folders within a single root directory which is created if necessary during the session setup (see Getting Started).

Under the root directory, there is a set of folders each corresponding to the data collected by a particular service. The service folder name is assigned during the particular service object creation and takes a default name if not specified (“Ensembl” for AnnotationEnsembl object “DAVID” for AnnotationDavid object etc.). It is possible to create multiple folders for a given service object type in case the folder name specified explicitly thus allowing to cache multiple instances of data from the

same online resource. This feature allows to store data available from the same service at different moments allowing to investigate how the given data source was evolving over time (see Use Case 2).

For a given service folder, there is a set of subfolders each of which corresponds to a particular micro array name (“HG-U133\_Plus\_2” etc.). Each micro array type related subfolder contains a set of files storing the R data objects (.RData files). There are two types of objects which are stored under the micro array folder. The first is a (single) data frame containing the raw data retrieved from a particular service, possibly containing multiple columns each for a particular data attribute. The second is a set of data frames representing the mapping between the particular attribute pairs (Uniprot Accession to probeset ID mapping as an example) extracted from a raw data frame during the last step of the ID mapping information retrieval process. The “two tier” caching allows to further speed up the data retrieval process for a particular attribute pair eliminating the need to read the entire raw data frame, whose size can easily exceed few dozen MB.

## 2.4 The output data formats

### 2.4.1 Identifier mapping data format

The identifier mapping related information (ID Map) is represented by a data frame consisting of two columns in character format. The first (primary) column contains a unique set of identifiers of a particular type (Uniprot accessions being an example) which represent the identifier mapping source. The second (secondary) column contains the results of mapping of a given (primary) identifier to the target identifier type (Affymetrix probeset ID as an example) obtained as a result of a query to the particular online service. As a single source identifier can map to multiple target identifiers, the secondary column elements are represented by comma separated strings.

The column names can be set during the service object initialization or can be deduced from the raw data if this information is not provided. In later case, the ID Map column names will depend on the particular service in use. If a further comparison of ID mapping performance for various online services is desirable, it is recommended to set the unified column names during the service objects initialization so to the nature of ID Map data (“Uniprot” and “Affy” for Uniprot accessions to Affymetrix probeset ID mapping as an example). In further writing, we will refer to the names of a first and second columns as primary and secondary keys correspondingly.

### 2.4.2 Complete raw data set format

The system always returns the complete raw data set as a data frame. As each particular online service uses it’s own data format, the data frame format will be different for each particular service and the query type. If the detailed information on the data frame content is desirable, we refer user to the documentation from a particular online service provider

## 3 Examples

### 3.1 Session initialization and setup

The data retrieval session always starts with the call to the package initialization function which tells the system which directory should be used for data caching. If specified directory does not exist yet, it will be created during the initialization process. As some online services require credentials (user ID and password) to grant user an access to the data, the next step of initialization process for a newly created caching directory should be setting the credentials for a service which requires them. At the moment, the only service directly requiring credentials is an Affymetrix query-based service (AnnotationNetAffx), however the rest of implemented services also utilizes some of the AnnotationAffx functionality (micro array name and/or probeset ID list retrieval) and therefore requires setting the same credentials during initialization process. Note that setting credentials is mandatory only for newly created caching directory and can be omitted in subsequent data retrieval sessions. The final step in initialization process is setting the verbosity level, defining the micro array type, mnemonic names for source and target identifier types (primary and secondary key values) as well as location of

files and/or folders for file-based services. A typical initialization script for a data retrieval session is given below:

```
> library(IdMappingRetrieval)
> # setup verbosity level,array type and DB keys
> arrayType <- "HG-U133_Plus_2";
> primaryKey <- "Uniprot";
> secondaryKey <- "Affy";
> #initialize caching directory structure
> Annotation$init(directory="./annotationData",verbose=TRUE);

Creating annotation archive directory: ./annotationData ...

> #set Affymetrix credentials
> AnnotationAffx$setCredentials(
+   user="alex.lisovich@gmail.com",password="125438",verbose=TRUE);

Creating service archive directory: ./annotationData/Affymetrix ...
```

### 3.2 Service object initialization

A few examples of different service objects initialization are given below.

```
> #create Affymetrix service object
> annAffx_Q <- AnnotationAffx(cacheFolderName="Affymetrix",
+                             primaryColumn="Probe.Set.ID",
+                             secondaryColumn="SwissProt",
+                             swap=TRUE);
> #create Ensembl file-based service object
> annEnsemblCsv <- AnnotationEnsemblCsv(cacheFolderName="EnsemblCsv",
+                                       primaryColumn=c("UniProt.SwissProt.Accession",
+                                                         "UniProt.TrEMBL.Accession"),
+                                       secondaryColumn=NA,
+                                       swap=FALSE, full.merge=TRUE,
+                                       df_filename="ENSEMBL_biomart_export.txt");
> #create Envision service object
> annEnvision <- AnnotationEnvision(cacheFolderName="EnVision",
+                                   primaryColumn=c("UniProt.SwissProt.Accession",
+                                                     "UniProt.TrEMBL.Accession"),
+                                   secondaryColumn=NA,
+                                   swap=TRUE, species="Homo sapiens", full.merge=TRUE);
```

### 3.3 Collecting ID maps from multiple services

Typically, the IdMappingRetrieval package is used in conjunction with the IdMappingAnalysis package which accepts the set of collected identifier mapping data in the form of IdMap object list . The steps involved in the process include:

- Defining array type and primary/secondary key values for a resulting IdMap object set \*.
- Defining the location of a file or directory containing preloaded data for file-based services
- Creating a set of service retrieval objects, one per each particular service
- Retrieving the Id Map through the getIdMap() method or a complete raw data frame through the getDataFrame() method.

\*The need to define the primary/secondary column names for resulting IdMap object set on top of primaryColumn/secondaryColumn arguments used during the service object creation arises from the fact that there is no standard naming convention for the raw data column names among the services. Therefore for a subsequent analysis of identifier mapping performance among the services we provide the arguments primaryKey and secondaryKey to getIdMap() and getIdMapList() methods for Annotation and ServiceManager objects respectively.

The script below illustrates the process of collecting such data from all services implemented within the IdMappingRetrieval package.

```
> # create service objects
> services <- list();
> services$NetAffx_F <- AnnotationAffx("Affymetrix");
> services$DAVID_Q <- AnnotationDavid("DAVID",species="Homo sapiens");
> #collect Id Map data
> idMapList <- lapply(services, getIdMap,
+                     arrayType=arrayType,
+                     primaryKey=primaryKey,
+                     secondaryKey=secondaryKey,
+                     force=FALSE, verbose=TRUE);
```

retrieving master annotation table...

Creating directory for array specific RData...

Reading compressed csv file...

Saving annotation data...

Swapping ID Map keys

unfolding ID Map ...

converting unique pairs into the ID Map ...

Saving ID map

Creating service archive directory: ./annotationData/DAVID ...

retrieving master annotation table...

Creating directory for array specific RData...

Retrieving Affymetrix Probeset List for HG-U133\_Plus\_2 ...

retrieving master annotation table...

Loading annotation data...

Saving ID map

Retrieving ID conversion choices...

URL= <http://david.abcc.ncifcrf.gov/tools.jsp> :

Uploading the ID list for AFFYMETRIX\_3PRIME\_IVT\_ID

URL= <http://david.abcc.ncifcrf.gov/conversion2.jsp> :

Request conversion to UNIPROT\_ACCESSION

Retrieving conversion file: [http://david.abcc.ncifcrf.gov/data/download/conv\\_D683E8366980142923802722](http://david.abcc.ncifcrf.gov/data/download/conv_D683E8366980142923802722)

Saving annotation data...

converting unique pairs into the ID Map ...

Swapping ID Map keys

unfolding ID Map ...

converting unique pairs into the ID Map ...

Saving ID map

>

This initial retrieval takes some time. But repeating the retrieval is fast due to caching.

```

> idMapList <- lapply(services, getIdMap,
+                     arrayType=arrayType,
+                     primaryKey=primaryKey,
+                     secondaryKey=secondaryKey,
+                     force=FALSE, verbose=TRUE);

Loading ID map...
Loading ID map...

> names(idMapList);

[1] "NetAffx_F" "DAVID_Q"

> idMapList$NetAffx_F[100:110,];

      Uniprot                                     Affy
F5GWT5  F5GWT5                                1487_at,203193_at
F5H0E9  F5H0E9                                1487_at,203193_at
H0YGT3  H0YGT3                                1487_at,203193_at
Q569H8  Q569H8                                1487_at,203193_at
Q6GMQ8  Q6GMQ8                                1487_at,203193_at
Q8N4S8  Q8N4S8                                1487_at,203193_at
Q96F89  Q96F89                                1487_at,203193_at
Q96IO2  Q96IO2                                1487_at,203193_at
P11509  P11509 1494_f_at,207244_x_at,207718_x_at,211295_x_at,214320_x_at
MOR2Z4  MOR2Z4          1494_f_at,207244_x_at,211295_x_at,214320_x_at
P20853  P20853          1494_f_at,207718_x_at

> idMapList$DAVID_Q[100:110,];

      Uniprot                                     Affy
O75488  O75488                                1569927_at,209421_at
Q6VBB7  Q6VBB7                                1569927_at,209421_at
Q0ZAI9  Q0ZAI9                                1569927_at,209421_at
O14754  O14754 203248_at,1554045_at,212534_at,203247_s_at,242210_at
Q6ICR5  Q6ICR5 203248_at,1554045_at,212534_at,203247_s_at,242210_at
Q53YE4  Q53YE4 203248_at,1554045_at,212534_at,203247_s_at,242210_at
P17028  P17028 203248_at,1554045_at,212534_at,203247_s_at,242210_at
Q8IZN4  Q8IZN4 203248_at,1554045_at,212534_at,203247_s_at,242210_at
Q9NR58  Q9NR58 217018_at,219307_at,232788_at,1555886_at,216669_at
Q5VU38  Q5VU38 217018_at,219307_at,232788_at,1555886_at,216669_at
Q5VU39  Q5VU39 217018_at,219307_at,232788_at,1555886_at,216669_at

```

Instead of map objects, the results can be stored as data frames.

```

> #collect complete data frames
> dfList <- lapply(services, getDataFrame,
+                 arrayType=arrayType,
+                 force=FALSE, verbose=TRUE);

Loading annotation data...
Loading annotation data...

> names(dfList);

[1] "NetAffx_F" "DAVID_Q"

```

```
> names(dfList$NetAffx_F);

[1] "Probe.Set.ID"                "GeneChip.Array"
[3] "Species.Scientific.Name"     "Annotation.Date"
[5] "Sequence.Type"              "Sequence.Source"
[7] "Transcript.ID.Array.Design." "Target.Description"
[9] "Representative.Public.ID"    "Archival.UniGene.Cluster"
[11] "UniGene.ID"                 "Genome.Version"
[13] "Alignments"                 "Gene.Title"
[15] "Gene.Symbol"                "Chromosomal.Location"
[17] "Unigene.Cluster.Type"       "Ensembl"
[19] "Entrez.Gene"                "SwissProt"
[21] "EC"                         "OMIM"
[23] "RefSeq.Protein.ID"          "RefSeq.Transcript.ID"
[25] "FlyBase"                    "AGI"
[27] "WormBase"                   "MGI.Name"
[29] "RGD.Name"                   "SGD.accession.number"
[31] "Gene.Ontology.Biological.Process" "Gene.Ontology.Cellular.Component"
[33] "Gene.Ontology.Molecular.Function" "Pathway"
[35] "InterPro"                   "Trans.Membrane"
[37] "QTL"                        "Annotation.Description"
[39] "Annotation.Transcript.Cluster" "Transcript.Assignments"
[41] "Annotation.Notes"

> names(dfList$DAVID_Q);

[1] "From"      "To"      "Species"  "Gene.Name"
```

### 3.4 Collecting data from various services using service manager

An alternative way of achieving the same results is to use the *ServiceManager* object, a container for service objects. The *ServiceManager* class generates the full set of default services automatically. The following script repeats the work of the previous code, more compactly.

```
> #create service manager object containing set of default services
> svm <- ServiceManager(ServiceManager$getDefaultServices());
> names(getServices(svm));

[1] "Affx_Q"      "Affx_F"      "DAVID_Q"      "DAVID_F"      "Ensembl_Q"
[6] "Ensembl_F"   "EnVision_Q"

> idMapList <- getIdMapList(svm,
+                           arrayType=arrayType,
+                           selection=c("NetAffx_Q", "DAVID_Q"),
+                           primaryKey=primaryKey,
+                           secondaryKey=secondaryKey, verbose=TRUE);

retrieving data for service *** DAVID *** :
Loading ID map...

> #collect complete data frames
> dfList <- getDataFrameList(svm,
+                             arrayType=arrayType,
+                             selection=c("NetAffx_Q", "DAVID_Q"), verbose=TRUE);

retrieving data for service *** DAVID *** :
Loading annotation data...
```



### 3.5 Collecting data from various services interactively using service manager

If micro array type and/or the paths for file-based services are not defined (set to NULL), ServiceManager present user with dialog(s) allowing selecting it interactively. Note that in a code snippet below the set of default services is generated internally by invoking a static version of getIdMapList() and getDataFrameList().

```
> #collect complete IdMap data interactively
> #selecting array type and services to collect data from
> idMapList <- ServiceManager$getIdMapList(
+           primaryKey= primaryKey,
+           secondaryKey=secondaryKey, verbose=TRUE);
> #collect complete data frames interactively
> #selecting array type and services to collect data from
> dfList <- ServiceManager$getDataFrameList(
+           arrayType="menu",
+           selection="menu", verbose=TRUE);
```

## 4 Session information

This version of IdMappingRetrieval has been developed with R 2.11.0.

R session information:

```
> toLatex(sessionInfo())
```

- R version 3.2.0 RC (2015-04-08 r68161), x86\_64-apple-darwin10.8.0
- Locale: en\_US.UTF-8/en\_US.UTF-8/en\_US.UTF-8/C/en\_US.UTF-8/en\_US.UTF-8
- Base packages: base, datasets, graphics, grDevices, methods, stats, utils
- Other packages: bitops 1.0-6, IdMappingRetrieval 1.16.0, R.methodsS3 1.7.0, R.oo 1.19.0, rChoiceDialogs 1.0.6, RCurl 1.95-4.5, rJava 0.9-6, XML 3.98-1.1
- Loaded via a namespace (and not attached): AffyCompatible 1.28.0, AnnotationDbi 1.30.0, Biobase 2.28.0, BiocGenerics 0.14.0, biomaRt 2.24.0, Biostrings 2.36.0, DAVIDQuery 1.28.0, DBI 0.3.1, ENVISIONQuery 1.16.0, GenomeInfoDb 1.4.0, IRanges 2.2.0, parallel 3.2.0, RSQLite 1.0.0, S4Vectors 0.6.0, stats4 3.2.0, tools 3.2.0, XVector 0.8.0, zlibbioc 1.14.0

## 5 Acknowledgments

We would like to thank Rafael Jimenez from the Enfin EnCore team for an enjoyable discussion and the great effort he has put into expanding the ID mapping related EnCore web services to better suite our needs. We also would like to thank Himanshu Grover, a graduate student in Biomedical Informatics at the University of Pittsburgh, for the great effort and useful suggestions he has made during the package preparation and testing.

## 6 References

- [1] Liu, G., NetAffx: Affymetrix probesets and annotations. Nucleic Acids Res. 2003, 31(1):82-6.
- [2] Huang D.W., Sherman B.T., Tan Q., Kir J., Liu D., Bryant D., Guo Y., Stephens R., Baseler M.W., Lane H.C. et al. (2007) DAVID Bioinformatics Resources: expanded annotation database and novel algorithms to better extract biology from large gene lists. Nucleic Acids Res., 35, W169-W175.

- [3] Huang D.W., Sherman B.T. and Lempicki R.A. (2008) Systematic and integrative analysis of large gene lists using DAVID bioinformatics resources. *Nat. Protoc.*, doi: 10.1038/nprot.2008.211.
- [4] Xosé M. Fernández-Suárez and Michael K. Schuster Using the Ensembl Genome Server to Browse Genomic Sequence Data. UNIT 1.15 in *Current Protocols in Bioinformatics*, Jun 2010. [www.ncbi.nlm.nih.gov/pubmed/20521244](http://www.ncbi.nlm.nih.gov/pubmed/20521244)
- [5] Giulietta M Spudich and Xosé M Fernández-Suárez Touring Ensembl: A practical guide to genome browsing *BMC Genomics* 2010, 11:295 (11 May 2010)
- [6] Florian Reisinger, Manuel Corpas, John Hancock, Henning Hermjakob, Ewan Birney and Pascal Kahlem. ENFIN - An Integrative Structure for Systems Biology. *Data Integration in the Life Sciences Lecture Notes in Computer Science*, 2008, Volume 5109/2008, 132-143, DOI: 10.1007/978-3-540-69828-9\_13
- [7] Henrik Bengtsson, The R.oo package - Object-oriented programming with references using standard R code, Talk at the DSC 2003 conference, Vienna, March 22, 2003.

## Appendix A

### Obtaining the Affymetrix ID Mapping data files using the NetAffx batch query system and HG-U133\_Plus\_2 array as an example.

1. Create a set of the files for upload containing the set of Affymetrix Probe Set IDs for a given array ensuring the number of IDs for each file does not exceed 10000 (NetAffx batch query limit):

```
AnnotationNetAffx$createSubmission(
  chunk=10000,folder="/NetAffx_submission",verbose=TRUE);
```

2. Register on the Affymetrix Web site ([www.affymetrix.com](http://www.affymetrix.com)) if you do not have an account yet.
3. Select "NetAffx" item in the menu at the top of the page and then select "Batch Query" item in "3' IVT Expression" tab on the NetAffx page. The "Batch Query" page will popup.
4. Create a custom annotation view which should include Probe Set ID as primary (first) identifier and SwissProt as a secondary if not done yet as follows: - Select "Custom Annotation Views" in the leftmost panel of the page (NetAffx Analysis Center). The "Create/Edit Custom Views" page will popup. - Enter the custom view name ("Affy2Uniprot" for example) and select Probe Set ID and SwissProt items in "Select Fields for your view" list. Click Submit. - Go back to the "Batch Query" page. Select "Human Genome U133 Plus 2.0 array" and "Probe Set ID" in "Select a GeneChip Array" and "Select the search type" list boxes correspondingly.
5. Select the first upload file from the set created at step 1 using "Browse" button and click "Submit". The "Show Results" page will popup.
6. Click on "Export Results" item. The "Export Center" page will popup. Select "Tab Separated Value" and "Affy2Uniprot" in "Select File Format" and "Select a View" list boxes correspondingly. Click "Submit" and then save the exported file when prompted.
7. Repeat steps 3, 5 and 6 for files created at step 1. Copy the set of resulting .tsv files into directory of your choice and use them when working with AnnotationNetAffx service retrieval object.
8. After downloading the file, the ID Map data can be retrieved by creating the AnnotationNetAffx object with df\_filename slot containing the path to the downloaded file directory (or the character vector of paths to file set) and then calling Annotation.getIdMap() on this object. If df\_filename is not specified (NULL by default) the File Open dialog will popup during the data retrieval allowing to specify the file set interactively.

## Appendix B

**Obtaining the DAVID ID Mapping data files through the DAVID knowledge base request form using Homo Sapiens as species and AFFY\_ID and UNIPROT\_ID as primary and secondary keys.**

1. Using your Web Browser, open the following url:  
[http://david.abcc.ncifcrf.gov/knowledgebase/DAVID\\_knowledgebase.html](http://david.abcc.ncifcrf.gov/knowledgebase/DAVID_knowledgebase.html)
2. Click on the “Request data from the DAVID Knowledgebase!! “ web page item. The “DAVID KNOWLEDGEBASE LOGIN” web page will open. If you are not registered yet, follow the “registration page” link to register. Enter your email address used during the registration process and then click “Login”.
3. The license agreement page “NATIONAL INSTITUTES OF HEALTH ACADEMIC, GOVERNMENT, NON-PROFIT RIGHT TO USE SOFTWARE AGREEMENT for DAVID Knowledgebase 1.0” will open. Read the agreement and then click “Accept”.
4. The “DAVID Knowledgebase Request” page will open. Following the instructions define species as “Homo Sapiens” (Step 1, Select Species), AFFY\_ID as a primary key (Central Identifier, Step 2) and UNIPROT\_ID as secondary key (Select Annotation Categories/Main Accessions, Step 3). Click “Submit”.
5. The confirmation page containing the summary of your request will open. In couple days you will receive notification by email containing the link to the download page.
6. After downloading the file, the ID Map data can be retrieved by creating the AnnotationDavidCsv object with df\_filename slot containing the path to the downloaded file and then calling Annotation.getIdMap() on this object. If df\_filename is not specified (NULL by default) the File Open dialog will popup during the data retrieval allowing to specify the file location interactively.

## Appendix C

**Obtaining the Ensembl ID Mapping data file using the HG-U133\_Plus\_2 array as an example.**

1. Open the Ensembl main page ([www.ensembl.org](http://www.ensembl.org)) and then select “BioMart” item from the menu at the top.
2. Choose “Ensembl Genes 60” as a database and select “Homo sapiens genes (GRCh37.p2)” as a data set.
3. Define the attributes by clicking on the “Attributes” item in the left panel and then check the following items in the right panel: - “GENE/ Ensembl Gene ID”, - “External References/ UniProt/TrEMBL Accession” , - “External References/ UniProt/SwissProt Accession” and - “Microarray/ Affy HG U133-PLUS-2” Click on “Results” button above the left panel. The Results will be displayed in the right panel.
4. In the Results panel, set “Export all results to” as “File” and “CSV” and check the box “Unique Results only”. Press “GO” and then save the file when prompted.
5. After downloading the file, the ID Map data can be retrieved by creating the AnnotationEnsembleCsv object with df\_filename slot containing the path to the downloaded file and then calling Annotation.getIdMap() on this object. If df\_filename is not specified (NULL by default) the File Open dialog will pop up during the data retrieval allowing to specify the file location interactively.