# Package 'lumiBarnes'

March 28, 2024

Type Package	
Title Barnes Benchmark Illumina Tissues Titration Data	
Version 1.42.0	
<b>Date</b> 2010-10-12	
Author Pan Du	
Maintainer Pan Du <dupan@northwestern.edu></dupan@northwestern.edu>	
<b>Description</b> The Barnes benchmark dataset can be used to evaluate the algorithms for Illumina microarrays. It measured a titration series of two human tissues, blood and placenta, and includes six samples with the titration ratio of blood and placenta as 100:0, 95:5, 75:25, 50:50, 25:75 and 0:100. The samples were hybridized on Human-Ref-8 BeadChip (Illumina, Inc) in duplicate. The data is loaded as an LumiBatch Object (see documents in the lumi package).	
License LGPL	
<b>Depends</b> R (>= 2.0), Biobase (>= 2.5.5), lumi (>= 1.1.0)	
biocViews ExperimentData, Tissue, MicroarrayData, ChipOnChipData	
git_url https://git.bioconductor.org/packages/lumiBarnes	
git_branch RELEASE_3_18	
git_last_commit 8d90ec5	
git_last_commit_date 2023-10-24	
Repository Bioconductor 3.18	
Date/Publication 2024-03-28	
R topics documented:	
lumiBarnes	4
Index	•

2 lumiBarnes

lumiBarnes

Barnes Benchmark Illumina Tissues Titration Data

## Description

The Barnes data set measured a titration series of two human tissues, blood and placenta. It includes six samples with the titration ratio of blood and placenta as 100:0, 95:5, 75:25, 50:50, 25:75 and 0:100. The samples were hybridized on HumanRef-8 BeadChip (Illumina, Inc) in duplicate. See (Barnes, et al., 2005) for details. The data is saved as a LumiBatch object and should be use together with lumi package.

Because the Barnes data utilized the pre-released version of HumanRef-8 version 1 BeadChip, some probes on the chip do not exist in the public released HumanRef-8 version 1 BeadChip. For annotation consistence, these probes was removed in the lumiBarnes package. For the interested users, the raw data can be downloaded from the paper companion website: http://www.bioinformatics.ubc.ca/pavlidis/lab/platformCom

### Usage

data(lumiBarnes)

#### **Format**

lumiBarnes is a LumiBatch-class object.

## Source

Barnes, M., Freudenberg, J., Thompson, S., Aronow, B. and Pavlidis, P. (2005) Ex-perimental comparison and cross-validation of the Affymetrix and Illumina gene expression analysis platforms, Nucleic Acids Res, 33, 5914-5923.

## **Examples**

data(lumiBarnes)
lumiBarnes

## **Index**

 $*\ datasets$ 

lumiBarnes, 2

 ${\tt lumiBarnes}, {\tt 2}$